


# Migrated hernia Rake®

## 【migrated hernia rake for TF】

Migrated hernia rake for TF has one jaw that rotates in a plane perpendicular to the mouth opening. The far-migrated hernia is held between the sawtooth surface of the jaw mouth and the dorsal side of the vertebral body, while then it is pulled out. The herniated disc was removed in the PELD transforaminal approach as follows: the jaw tip was advanced toward the posterior side of the migrated hernia under the posterior longitudinal ligament or dura. The herniated disc was held from the posterior side with the jaw by rotating the shaft of the instrument on its axis, and it was pulled out or stretched to the proximal side by pulling the instrument while rotating the jaw back to the straight position, and subsequently removed using conventional forceps. This maneuver is safe even if in blind manipulation, because it does not bite the tissue, but pinches in a safe direction. The jaw is long (9 or 12 mm) and moves in a right or left direction.



migrated hernia rake® for TF	Type	Overhang Length	OD	Working Length
	R (Right)	9mm	3.5mm	360mm
	L (Left)			
	R (Right)	12mm		
	L (Left)			


※The photo is not actual size. The size is subject to change without notice.  
Various custom sizes are available.

Device name : フェジ鉗子  
 Registration Number : 21B2X10011000001  
 Device Class : Class I  
 Adachi Industry Co.,Ltd.  
 TEL : +81-575-22-0052  
 Vol.1-2022.10

# migrated hernia rake®

## 【migrated hernia rake for all-around】

Migrated hernia rake for all-around has one jaw that rotates in the direction of the mouth opening. The function is the same as that for TF. It can be used in transforaminal, posterolateral, and interlaminar approach.



migrated hernia rake® for all-around	Type	Overhang Length	OD	Working Length
	Up	10mm	3.5mm	360mm
	Down			

※The photo is not actual size. The size is subject to change without notice.  
Various custom sizes are available.

Device name : フェジ鉗子  
 Registration Number : 21B2X10011000001  
 Device Class : Class I  
 Adachi Industry Co.,Ltd.  
 TEL : +81-575-22-0052  
 Vol.1-2022.10