







# Working sheath for percutaneous full-endoscope



Wide pitch screw for the clearance between the outer sheath and the endoscope allows for quick and stepless adjustment.

Clearance adjustment screw (removed.)

	straight
	75° oblique
	35° ~45° oblique
	duckbill
	stepped

OD	ID	length
9.5mm	8.6mm	112mm
8.0mm	7.1mm	170mm
7.5mm	6.5mm	170mm

※ We can produce it to accommodate the size of full-endoscopes of various manufacturers. Please contact us for details.

※ The photo is not actual size. The size is subject to change without notice. Various custom sizes are available.

Device name : 脊椎全内視鏡用外筒  
 Registration Number : 21B2X10011000010  
 Device Class : Class I  
 Adachi Industry Co., Ltd.  
 TEL : +81-575-22-0052  
 Vol. 1-2022. 11 (P 1 / 1)